

# myTAP Vertical Shim

## Instructions



## INTENDED USE

The **myTAP** vertical shims are an optional comfort accessory to the **myTAP** oral appliance. They increase the vertical opening of the mouth and come in small, medium and large sizes that progressively increase the opening of the mouth.

## ASSEMBLY

To assemble remove the lower tray from the upper tray and slide the shim over the grey threaded shaft of the lower tray. The sizing letter (S, M or L) should be located on top of the lower tray orientated towards the mouth. Reassemble the upper and lower trays of your **myTAP**, back to the original position, before use.

Only one shim should be used at a time.

## CLEANING

Refer to **myTAP** instructions for cleaning and warnings.